

## Type of Conventional Child Safety Seat by Ambulance Installation Location

	Infant-only	Convertible	Forward-Facing Only	Combination Seat		Belt Positioning Booster	Dream Ride Car Bed	Other Specialized Restraints
				Harness Mode	Booster Mode			
<b>Front Passenger Seat*</b>	Yes‡‡ - CSS rear facing per CSS manufacturer instructions.	Yes - Forward or (rear facing‡‡) per child's weight & CSS mfr instructions.	Yes - per CSS manufacturer instructions.	Yes - per CSS manufacturer instructions.	Yes - per CSS manufacturer instructions. Lap/shoulder belt required.	Yes - per CSS manufacturer instructions. Lap/shoulder belt required.	No	Yes‡‡ - per manufacturer instructions. Tether strap required.
<b>Capt's Chair <u>Locked in Forward Facing Position</u>**‡</b>	Yes‡‡ - CSS rear facing per CSS manufacturer instructions.	Yes‡‡ - Forward or rear facing per child's weight and CSS manufacturer instructions.	Yes‡‡ - per CSS manufacturer instructions.	Yes‡‡ - per CSS manufacturer instructions.	Yes‡‡ - per CSS manufacturer instructions. Lap/shoulder belt required.	Yes‡‡ - per CSS manufacturer instructions. Lap/shoulder belt required.	No	No
<b>Capt's Chair <u>Locked in Rear Facing Position</u></b>	No	Yes - Face CSS toward rear of ambulance. <b>Use belt path marked "forward facing"</b> .	Yes -Face CSS toward rear of ambulance.	Yes -Face CSS toward rear of ambulance.	No	No	No	No
<b>Cot/Stretcher <u>Best location for pediatric patient</u></b>	No	Yes - Face CSS toward rear of ambulance. Secure using 2 belts & both belt paths.	No	No	No	No	Yes - Place head away from sidewall. Requires 2 sets of bed loops and 2 belts.	No
<b>Side Facing Bench/Chair</b>	No	No	No	No	No	No	No	No

\* Airbags off and vehicle seat pushed back.

\*\* Check with ambulance manufacturer to determine if captain's chair has been crash tested in forward facing position.

‡ Open space in front of seat must allow for excursion of head and extremities.

‡‡ **Caution** - These locations and CSS should only be used if other options are unavailable.